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Nastro inchiostrato : 739209 (ERC09 INK)

TABLE OF CONTENTS

Declaration of conformity	page 17
Symbols used	page 18
Safety rules	page 19
Features	page 20
Description of parts	page 21-22
Dimensions	page 22
Installations	page 23
Operating instructions	page 24
Maintenance	page 26-27-28
Troubleshooting	page 29
Data sheet	page 30

16

"CE" Declaration of Conformity

In accordance with standards ISO/IEC Guide 22 and EN 45014

Manufacturer's name: Custom Engineering S.p.A.

Manufacturer's address: Strada Berettine 2

Fontevivo (Parma)

Italy

Declares that the product:

Product name: Printer

Product type: Steriprint

Sold by: FARO S.p.A.

Via Faro, 15 Ornago (Mi)

Italy

Is in conformity with the following directives:

Low Voltage directive 73/23/CEE; 93/68/CEE

Electromagnetic Compatibility Directive 89/336/CEE; 92/31/CEE; 93/68/CEE

In accordance with the following standards:

EN 55022

Class B	characteristics of radio disturbance produced by information technology equipment	1994
EN 61000-3-2	Limits for the emission of harmonics in current. Input power $\leq 50W$	1995
EN 61000-3-3	Limits for voltage fluctuations and flickers	1995
EN 50082-1	Electromagnetic compatibility - General immunity requirements. Part 1: Residential, commercial and light industry environments.	1992
EN 61000-4-2	Electrostatic discharge immunity tests. 4KVcontact discharge, 8KV air discharge	1995
EN 61000-4-4	Electrical fast transient/burstimmunity tests. Signal lines 0.5KV Power lines AC 1KV	1995
ENV 50140	Radio-frequency irradiated electromagnetic fields - Immunity test. 3V/m, 80MHz-1000MHz, 80% 1KHz AM	1993
EN 60950 A1+A2	Information technology equipment including electrical office appliances. Safety	1992 1993

Limits and methods of measuring the

1994

17

SYMBOLS USED IN THIS MANUAL

N.B.

Gives important information or suggestions relative to the use of the printer

WARNING

The information marked with this symbol must be carefully heeded to safeguard against damaging the printer

DANGER

The information marked with this symbol must be carefully heeded to safeguard against injury to the operator

SAFETY RULES

Read and keep the following instructions.

- Observe all warnings and follow all instructions attached to the printer.
- Before cleaning the printer, disconnect the feed cable.
- Clean the printer with a damp cloth. Do not use liquid or spray products.
- Do not operate the printer near to water.
- Do not place the printer on unsteady surfaces. It could fall and get seriously damaged.
- Do not place the printer on soft surfaces or in poorly ventilated environments.
- Position the printer in such a way as to ensure that the cables connected to it will not be damaged.
- Use the type of electricity supply marked on the printer label. In the event of uncertainty, contact the seller.
- Ensure that the printer's electricity supply is grounded and that it is protected by a differential switch.
- If the printer is fed through an extension lead, ensure that the total absorption of the equipment connected to it does not exceed the maximum admissible current for that type of extension and that it does not, in any event, exceed 15 amperes.
- Do not obstruct the vents.
- Do not put objects of any kind inside the printer as they could cause a short circuit or damage parts which could affect its performance.
- Do not spill liquids on the printer.
- Do not carry out technical operations on the printer with the exception of the scheduled maintenance operations specifically indicated in the user's manual.

FEATURES

The STERIPRINT is a practical, user-friendly desk printer.


It is, therefore, the ideal solution for applications which require the immediate printing of data, whether they be of an industrial, professional or laboratory nature. Typical fields of application are: POS, weighing systems, receipts (not for tax purposes) as well as for security, controlling and diagnostics purposes. It is equipped with a rapid, 8-needle impact print mechanism and uses 57.5 mm paper rolls.

The printer is distinguished by the fact that it has an internal power supply and is extremely compact.

It has an RS232 serial interface. It is equipped with a Real Time Clock.


DESCRIPTION OF PARTS

- 1- Cover
- 2- Paper outfeed
- 3- Feed key
- 4- Printer base
- 5- Printer body
- 6-LED
- 7- Print key
- 8- Print mechanism
- 9- Keypad


DESCRIPTION OF PARTS

- 1- Serial connector
- 2- ON/OFF switch
- 3- Feed connector

DIMENSIONS


INSTALLATION

Remove the printer from the box, taking care not to damage the packing material, as it could be needed for future transportation of the machine.

Ensure that all the components illustrated are in fact present and that they are in perfect condition. If this is not the case, contact the after-sales assistance department immediately.

- 1. Upper tray
- 2. Instruction Manual
- 3. Feed cable
- 4. Suspended upper packing
- 5 Printer
- 6. Paper roll (inside the printer)
- 7. Suspended lower packing
- 8. Lower tray
- 9. Case
- 10 Connection cable

- Unpack the printer
- Remove the protective tray
- Remove the upper suspended packing and remove the feed cable and manual
- Take the printer out of its bag.
- Keep the box, trays and suspended packing materials in the event of the having to be shipped to another destination.

The STERIPRINT has an RS232 serial interface with a 9-pin female connector. For the arrangement of the signals on the connector pins and for hooking up to the PC, please refer to the following table:

The printer has a feed cable 1.5 m long with a standard plug.

WARNING:

ensure that the mains voltage is compatible with the characteristics of the equipment.

- 1) Before connecting the printer to the mains, ensure that the power supply or system ON/OFF switch is in the OFF position.
- 2) Do not print without paper or ink ribbon; this leads to rapid deterioration of the needles.
- 3) Do not pull the printer carriage manually.

OPERATING INSTRUCTIONS

- Press the main switch (fig. 9 page 22).
- The green led lits on (when is on it indicates that the printer is working well).
- -The printer is ready to receive data from the autoclave.
- -FARO
- -Ornago (Mi)
- -SN..ST....
- -Program:
- -Cycle nr
- -Start:
- -(ora data)

Press	Temp	Time
-kPa	°C	min
-Fract. Vacuum		
-Ster.		
-Discarge		

-End:

-Dry

(ora - data)

-Press Balance

- OK

Function keys

PRINT: Button for programming

FEED key. When this is pressed the paper feeds forward manually. If this key is pressed briefly, when the RTCK option is installed, the date and time of day is printed.

Adjusting the clock

The time and date can be adjusted using the PRINT and FEED keys on the printer's front panel. To set, proceed as follows:


- While holdind down the FEED key, press the PRINT Key. The printer will print the time and date with an arrow indicating the digit to modified;
- Each time the PRINT key is pressed, the digit marked by the arrow will increase and an updated version will be printed.
- To proceed to modify another digit, press the FEED key again. Each time the printer will print the updated time and date, highlighting with an arrow the currently selected digit.
- To terminate the setting procedure, press PRINT and FEED at the same time, or scroll all the parameters

MAINTENANCE

Changing the paper roll

To change the paper roll, proceed as follows:


- 1) Open the upper cover and position the paper roll so that it rotates in the right direction, as shown in the figure;
- 2) Insert the end of the roll in the slit of the print mechanism;
- 3) Press the FEED key; a few centimetres of paper automatically feed out of the printer;
- 4) Tear off the paper and re-close the cover.

Replacing the ink ribbon

To replace the ink ribbon, proceed as follows:

- 1) Open the door and remove the used cartridge by pressing at the point marked PUSH (1), as shown in the figure;
- 2) Fit the new ribbon, ensuring that it is correctly positioned;
- **3)** Tighten the ribbon by rotating the knurled knob (2) in the direction indicated by the arrow; then re-close the door.

Paper jam

- 1) To remove the paper from the print mechanism, pull it upwards;
- 2) If the paper comes out the front, obstructing the printing area, **first** tear off the excess paper then carefully remove the scrap which has jammed.


Cleaning

To clean the printer, use a vacuum cleaner or a soft cloth.

Before cleaning the printer, disconnect the feed cable from the mains socket.

Do not use alcohol, solvents or hard-bristled brushes.

Do not allow water or other liquids to wet the internal mechanisms of the printer.

TROUBLESHOOTING

Disconnect the printer from the electricity supply and have it repaired by a specialized technician should any of the following conditions occur:

- The feed connector has been damaged.
- LIQUID has penetrated to the inside of the printer;
- The printer has been exposed to rain or water;
- The printer is not operating normally despite the instructions in the user's manual having been followed.
- The printer has been dropped and its case damaged.
- The performance of the printer is poor.
- The printer does not work.

DATA SHEET

Power supply : 230V ac 50 Hz

Environmentals conditions

Operating temperature : $0^{\circ} + 50^{\circ}$ C Operating humidity : -25° C $+70^{\circ}$ C Storage temperature : 10% - 90%

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GUARANTEE CERTIFICATE

GB

FARO grants to the final customer a 12-months guarantee starting from the date of the purchase invoice. Repairs under the guarantee shall be done at the FARO premises or by an authorized FARO Service person with transportation costs and risks at the expense of the purchaser. The guarantee certificate found on the endleaf of this booklet is considered valid only if it is entirely completed and bears the reseller's stamp. Repairs under guarantee (during the period of effectiveness) will be done only if the article to be repaired is accompanied by the shipping note or purchase invoice. The guarantee covers failures due to poor material quality or manufacturing detects and in case of legitimate claim the guarantee will allow repair or replacement free of charge. No compensation for damages and/or interest will be recognized. The guarantee is not deemed valid in the final judgement of FARO in case of tampering. damage, incorrect use, faulty maintenance or normal wear. The guarantee does not cover bulbs and fuses

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FARO S.p.A. reserves the right to change the specifications of this equipment without notice.

Steriprint

MANUALE DI ISTRUZIONI INSTRUCTION MANUAL MODE D'EMPLOI GEBRAUCHS ANWEISUNGEN MANUAL DE MANEJO

FARO S.p.A. Via Faro, 15

20060 Ornago (MI) Italia Tel. : 039 / 6878.1 (r.a.) Fax : 039 / 6010540 Internet: http://www.faro.it

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FARO S.p.A. si riserva il diritto di modificare, senza preavviso, le caratteristiche indicate nel presente manuale.